

Instruction manual



- 1.Name: Dot nevus scanning pen
- 2.Material: ABS
- 3.Input voltage: 5V=1A
- 4.Battery capacity: 1200mAh
- 5.Machine color:White, gold, silver, pink
- 6.Press and hold the start button for 3-5 seconds to turn on, and press and hold for 3-5 seconds to turn off.
- 7.When starting the machine, click on upshift (a total of 9 gears), double-click to turn on and off lighting.
- 8."ON"Click: Work once, long press to continue working.

Features: Portable and Rechargeable

How to use



Choose the right needles
Long press for two seconds
to open



Adjust the intensity



Working

Product include:

1* Dot nevus scanning pen

1* USB cable

1* manual

1* box

Manual de instrucciones



1. Nombre: lápiz de escaneo Dot Nevus
2. Material: ABS
3. Voltaje de entrada: 5V = IA
4. Capacidad de la batería: 1200mAh
5. Color de la máquina: blanco, dorado, plateado, rosa.
6. Mantenga presionado el botón de inicio durante 3-5 segundos para encender, y presione y mantenga presionado durante 3-5 segundos para apagar.

8."ON"Haga clic: Trabajar una vez, mantenga presionado para continuar trabajando

Características: portátil y recargable

Cómo utilizar



Elige las agujas adecuadas



Mantenga pulsado durante dos segundos
abrir



Ajustar la intensidad



Laboral

El producto incluye:

1 * bolígrafo de escaneo Dot nevus